

Anti-Pgm1 Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - P36871

Item Number: CM0003322



Introduction

Rabbit monoclonal antibody against human PGM1, validated for WB, IF/ICC, and ELISA. Reacts with human, mouse, and rat PGM1; observed molecular weight ~61kDa. Ideal for glucose metabolism and glycogen storage disease research.

[Learn More](#)

Parameter	Value
Product Name	PGM1 Rabbit mAb
Remarks/Alias	CDG1T; GSD14; PGM1
Species	Human
Gene ID (Human)	5236
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 465-562 of human PGM1 (NP_002624.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P36871
Protein Weight	61kDa
Shipping	Ice bag