

Anti-Phd2/Egln1 Rabbit Monoclonal Antibody For Wb, Elisa - Q9Gzt9

Item Number: CM0002656



Introduction

Anti-PHD2/EGLN1 rabbit monoclonal antibody, validated for Western blot (WB) and ELISA. Detects human PHD2/EGLN1 (HIF-PH2), a 46 kDa oxygen sensor involved in HIF- α hydroxylation and hypoxia signaling. Suitable for research use.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] PHD2/EGLN1 Rabbit mAb
Remarks/Alias	HPH2; PHD2; SM20; ECT3; HALAH; HPH-2; HIFPH2; ZMYND6; C1orf12; HIF-PH2; N1
Species	Human
Gene ID (Human)	54583
Gene ID	54583
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human PHD2/EGLN1 (NP_071334.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9GZT9
Protein Weight	46kDa
Shipping	Ice bag