

Anti-Phd3 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9H6Z9

Item Number: CM0030251



Introduction

Rabbit polyclonal antibody targeting human PHD3/EGLN3 (Q9H6Z9). Validated for WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat. Ideal for studying hypoxia, apoptosis, ferroptosis, and DNA damage response.

[Learn More](#)

Parameter	Value
Product Name	PHD3 Rabbit pAb
Remarks/Alias	PHD3; HIFPH3; HIFP4H3
Species	Human
Gene ID (Human)	112399
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human PHD3 (NP_071356.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H6Z9
Protein Weight	27kDa
Shipping	Ice bag