

Anti-Phf11 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9Uil8

Item Number: CM0026311

Introduction

PHF11 rabbit polyclonal antibody for detection of PHD finger protein 11, a Th1 cytokine regulator. Validated in WB, IHC-P, IF/ICC, ELISA; cross-reacts with human, mouse, rat.

[Learn More](#)

Parameter	Value
Product Name	PHF11 Rabbit pAb
Remarks/Alias	APY; BCAP; IGEL; IGER; IGHER; NYREN34; NY-REN-34; PHF11
Species	Human
Gene ID (Human)	51131
Gene ID	51131
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-292 of human PHF11 (NP_001035533.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UIL8
Protein Weight	38kDa
Shipping	Ice bag