


Anti-Phkg2 Rabbit Polyclonal Antibody For Wb, Elisa - P15735

Item Number: CM0014508



Introduction

Rabbit polyclonal antibody targeting human PHKG2, a catalytic subunit of phosphorylase b kinase involved in glycogenolysis. Suitable for WB and ELISA; cross-reacts with human, mouse, and rat. Expected band ~46 kDa.

[Learn More](#)

Product Name	PHKG2 Rabbit pAb
Remarks/Alias	GSD9C; PHKG2
Species	Human
GeneID (human)	5261
GeneID	5261
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P15735
Protein Weight	46kDa
Shipping	Ice bag