

Anti-Phlda2 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q53Ga4

Item Number: CM0030346

Introduction

Rabbit polyclonal antibody against PHLDA2 (IPL/TSSC3), validated for WB, IF/ICC, and ELISA in human and mouse. Suitable for placenta growth, imprinted gene, and cancer research studies.

[Learn More](#)

Parameter	Value
Product Name	PHLDA2 Rabbit pAb
Remarks/Alias	IPL; BRW1C; BWR1C; HLDA2; TSSC3; PHLDA2
Species	Human, Mouse
GeneID (Human)	7262
GeneID	7262
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-152 of human PHLDA2 (NP_003302.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q53GA4
Protein Weight	17kDa
Shipping	Ice bag