

Anti-Phospho-Akt-S473 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P31749 / P31751 / Q9Y243

Item Number: CM0001927



Introduction

Phospho-specific rabbit monoclonal antibody against AKT1/AKT2/AKT3 phosphorylated at serine 473, validated for WB, IHC-P, and ELISA in human, mouse, and rat samples. Ideal for apoptosis, metabolism, and cancer research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Akt-S473 Rabbit mAb
Remarks/Alias	AKT1/AKT2/AKT3; Phospho-Akt-S473
Species	Human
GeneID (human)	207/208/10000
GeneID	207/208/10000
Immunogen	Synthetic peptide; A synthetic phosphorylated peptide around S473 of human AKT1 (NP_005154.2)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P31749/P31751/Q9Y243
Protein Weight	48kDa/55kDa/51kDa/54kDa
Shipping	Ice bag