

Anti-Phospho-Ampka1/Ampka2-T183/T172 Rabbit Polyclonal Antibody For Wb, Elisa - Q13131 / P54646

Item Number: CM0001493



Introduction

Phospho-specific rabbit polyclonal antibody against AMPK α 1/ α 2 phosphorylated at Thr183/Thr172. Suitable for WB and ELISA. Detects human, mouse, and rat AMPK. Ideal for energy metabolism research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-AMPK α 1/AMPK α 2-T183/T172 Rabbit pAb
Remarks/Alias	AMPK α 1/AMPK α 2; Phospho-AMPK α 1/AMPK α 2-T183/T172
Species	Human
GeneID (Human)	5562/5563
GeneID	5562/5563
Immunogen	Synthetic peptide
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13131/P54646
Protein Weight	64kDa/65kDa/62kDa
Shipping	Ice bag