


Anti-Phospho-Aurka-T288/Aurkb-T232/Aurkc-T198 Rabbit Polyclonal Antibody For Wb, Elisa - O14965 / Q96Gd4 / Q9Uqb9

Item Number: CM0000808



Introduction

Phospho-specific rabbit polyclonal antibody recognizing AURKA (T288), AURKB (T232), and AURKC (T198). Validated for Western blotting and ELISA; reacts with human and rat samples. Ideal for cell cycle and cancer research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-AURKA-T288/AURKB-T232/AURKC-T198 Rabbit pAb
Remarks/Alias	AURKA/AURKB/AURKC; Phospho-AURKA-T288/AURKB-T232/AURKC-T198
Species	Human
GeneID (Human)	6790/9212/6795
GeneID	6790/9212/6795
Immunogen	Synthetic phosphorylated peptide around T288 of human AURKA (NP_003591.2)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	O14965/Q96GD4/Q9UQB9
Shipping	Ice bag