


Anti-Phospho-Cdk2-T160 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P24941

Item Number: CM0001429



Introduction

Phospho-CDK2-T160 Rabbit Polyclonal Antibody (pAb) validated for WB, IHC-P, ELISA. Detects human, mouse, rat CDK2 phosphorylated at Thr-160. Ideal for cell cycle, DNA damage, and oncology research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-CDK2-T160 Rabbit pAb
Remarks/Alias	CDKN2; p33(CDK2); Phospho-CDK2-T160
Species	Human
Gene ID (Human)	1017
Gene ID	1017
Immunogen	Synthetic peptide A synthetic phosphorylated peptide around T160 of human CDK2 (NP_001789.2?/NP_001249.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P24941
Protein Weight	30 kDa/34 kDa
Shipping	Ice bag