

Anti-Phospho-Egfr-Y1173 Rabbit Polyclonal Antibody For Wb, Elisa - P00533

Item Number: CM0001089



Introduction

Phospho-EGFR-Y1173 rabbit polyclonal antibody validated for Western blot and ELISA. Detects human EGFR phosphorylated at tyrosine 1173, a key receptor tyrosine kinase involved in cell signaling and cancer.

[Learn More](#)

Parameter	Value
Product Name	Phospho-EGFR-Y1173 Rabbit pAb
Remarks/Alias	ERBB; ERBP; HER1; mENA; ERBB1; PIG61; NISBD2; Phospho-EGFR-Y1173
Species	Human
GeneID (Human)	1956
GeneID	1956
Immunogen	Synthetic Peptide: A phospho specific peptide corresponding to residues surrounding Y1173 of human EGFR
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P00533
Protein Weight	134kDa
Shipping	Ice bag