


Anti-Phospho-Gys1 (S641) Rabbit Monoclonal Antibody For Wb, Elisa - P13807

Item Number: CM0001910



Introduction

Phospho-specific rabbit monoclonal antibody targeting human GYS1 at phosphorylation site S641. Validated for WB and ELISA, cross-reacts with human and mouse. Suitable for glycogen metabolism and diabetes research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-GYS1-S641 Rabbit mAb
Remarks/Alias	GSY; GYS; Phospho-GYS1-S641
Species	Human
GeneID (human)	2997
GeneID	2997
Immunogen	A synthetic phosphorylated peptide around S641 of human Glycogen synthase 1 (GYS1) (P13807)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P13807
Protein Weight	84kDa
Shipping	Ice bag