


Anti-Phospho-Hdac2 (S394) Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q92769

Item Number: CM0001190



Introduction

Phospho-specific rabbit polyclonal antibody targeting HDAC2 at Ser394, validated for WB, IHC-P, and ELISA applications. Reacts with human, mouse, and rat samples, ideal for epigenetic and transcriptional regulation studies.

[Learn More](#)

Parameter	Value
Product Name	Phospho-HDAC2-S394 Rabbit pAb
Remarks/Alias	HD2; RPD3; YAF1; KDAC2; Phospho-HDAC2-S394
Species	Human
Gene ID (Human)	3066
Gene ID	3066
Immunogen	Synthetic peptide; a phospho-specific peptide corresponding to residues surrounding S394 of human HDAC2
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q92769
Protein Weight	55kDa
Shipping	Ice bag