


Anti-Phospho-Hsf1-S307 Rabbit Polyclonal Antibody For Wb, Elisa - Q00613

Item Number: CM0000729



Introduction

Phospho-specific rabbit polyclonal antibody targeting HSF1 at Ser307. Validated for WB and ELISA, reacts with human HSF1. Ideal for studying heat shock response and stress signaling pathways.

[Learn More](#)

Parameter	Value
Product Name	Phospho-HSF1-S307 Rabbit pAb
Remarks/Alias	HSTF1; Phospho-HSF1-S307
Species	Human
Gene ID (Human)	3297
Gene ID	3297
Immunogen	Synthetic phosphorylated peptide around S307 of human HSF1.
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q00613
Protein Weight	57kDa
Shipping	Ice bag