


Anti-Phospho-Hsp27/Hspb1-S78 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P04792

Item Number: CM0000501



Introduction

Phospho-HSP27/HSPB1-S78 Rabbit Polyclonal Antibody validated for WB, IF/ICC, ELISA. Detects human HSPB1 phosphorylated at Ser78. Ideal for stress response, chaperone function, and actin organization research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-HSP27/HSPB1-S78 Rabbit pAb
Remarks/Alias	CMT2F; HMN2B; HSP27; HSP28; Hsp25; SRP27; HS.76067; HEL-S-102; Phospho-HSP27/HSPB1-S78
Species	Human
GeneID (Human)	3315
GeneID	3315
Immunogen	Synthetic peptide
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P04792
Protein Weight	23kDa
Shipping	Ice bag