


Anti-Phospho-Map2K4 (S257/T261) Rabbit Polyclonal Antibody For Wb And Elisa - P45985

Item Number: CM0001436



Introduction

Phospho-specific rabbit polyclonal antibody recognizing MAP2K4 phosphorylated at Ser257 and Thr261. Validated for WB and ELISA in human, mouse, and rat samples. Ideal for stress signaling and apoptosis research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-MAP2K4-S257/T261 Rabbit pAb
Remarks/Alias	JNKK; MEK4; MKK4; SEK1; SKK1; JNKK1; SERK1; MAPKK4; PRKMK4; SAPKK1; SAPKK-1; Phospho-MAP2K4-S257/T261
Species	Human
Gene ID	6416
Immunogen	Synthetic peptide (A synthetic phosphorylated peptide around S257 & T261 of human MAP2K4 (NP_003001.1))
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P45985
Protein Weight	44kDa
Shipping	Ice bag