


Anti-Phospho-Nf-Kb P65/Rela-S536 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q04206

Item Number: CM0001723



Introduction

High-quality Phospho-NF-kB p65/RelA-S536 rabbit polyclonal antibody validated for Western blot, IF/ICC, and ELISA. Reacts with human, mouse, rat. Ideal for NF-kB signaling pathway research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-NF-kB p65/RelA-S536 Rabbit pAb
Remarks/Alias	p65; CMCU; NFKB3; AIF3BL3; Phospho-NF-kB p65/RelA-S536
Species	Human
GeneID	5970
Immunogen	Synthetic phosphorylated peptide around S536 of human RelA (NP_068810.3)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q04206
Protein Weight	60kDa
Shipping	Ice bag