


## Anti-Phospho-Nf-Kb P65/Rela (Ser468) Rabbit Monoclonal Antibody For Wb, Ip, Elisa - Q04206

Item Number: CM0000695



### Introduction

Rabbit monoclonal antibody targeting phosphorylated NF-kB p65/RelA at Ser468. Validated for WB, IP, ELISA in human, mouse, rat. Ideal for studying inflammation, immunity, and cancer signaling pathways.

[Learn More](#)

| Parameter        | Value                                                                      |
|------------------|----------------------------------------------------------------------------|
| Product Name     | Phospho-NF-kB p65/RelA-S468 Rabbit mAb                                     |
| Remarks/Alias    | p65; CMCU; NFKB3; AIF3BL3; Phospho-NF-kB p65/RelA-S468                     |
| Species          | Human                                                                      |
| GeneID (Human)   | 5970                                                                       |
| Immunogen        | A synthetic phosphorylated peptide around S468 of human RELA (NP_068810.3) |
| Source           | Rabbit                                                                     |
| Category         | Phosphorylated Antibodies                                                  |
| Application      | WB, IP, ELISA                                                              |
| Cross Reactivity | Human, Mouse, Rat                                                          |
| SWISS            | Q04206                                                                     |
| Protein Weight   | 60kDa                                                                      |
| Shipping         | Ice bag                                                                    |