


Anti-Phospho-Nrf2-S40 Rabbit Mab For Wb, If/Icc, Ihc-P, Elisa - Q16236

Item Number: CM0000880



Introduction

Phospho-NRF2-S40 Rabbit mAb detects phosphorylated NRF2 (Ser40) in human, mouse, rat. Suitable for WB, IF/ICC, IHC-P, ELISA. Ideal for oxidative stress and cytoprotection research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-NRF2-S40 Rabbit mAb
Remarks/Alias	NRF2; HEBP1; Nrf-2; IMDDHH; Phospho-NRF2-S40
Species	Human
Gene ID	4780 (Human)
Immunogen	Synthetic Peptide A synthetic phosphorylated peptide around S40 of human NRF2 (Q16236)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16236
Protein Weight	68 kDa
Shipping	Ice bag