


Anti-Phospho-Pdgfrb (Y1021) Rabbit Polyclonal Antibody For Wb, Elisa - P09619

Item Number: CM0000765



Introduction

Rabbit polyclonal antibody detecting phosphorylated PDGFRB at Y1021, validated for Western blot and ELISA in human, mouse, and rat samples. Ideal for studying PDGF signaling in vascular biology and cancer.

[Learn More](#)

Parameter	Value
Product Name	Phospho-PDGFRB-Y1021 Rabbit pAb
Remarks/Alias	IMF1; KOGS; IBGC4; JTK12; PDGFR; PENTT; CD140B; PDGFR1; PDGFR-1; Phospho-PDGFRB-Y1021
Species	Human
Gene ID (Human)	5159
Gene ID	5159
Immunogen	Synthetic Peptide: A synthetic phosphorylated peptide around Y1021 of human Phospho-PDGFRB-Y1021 (NP_002600.1).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P09619
Protein Weight	124kDa
Shipping	Ice bag