


Anti-Phospho-Polr2A Ctd-S7 Rabbit Pab For Wb, Ihc-P, Elisa, Chip - P24928

Item Number: CM0000708



Introduction

Phospho-specific Rabbit Polyclonal Antibody targeting the CTD Ser-7 phosphorylation of POLR2A. Validated for WB, IHC-P, ELISA, CHIP in human, mouse, rat. Essential for studying Pol II transcription regulation and RNA processing.

[Learn More](#)

Parameter	Value
Product Name	Phospho-POLR2A CTD-S7 Rabbit pAb
Remarks/Alias	RPB1; RPO2; POLR2; POLRA; RPBh1; RPOL2; NEDHIB; RplILS; hsRPB1; hRPB220; Phospho-POLR2A CTD-S7
Species	Human
GeneID (Human)	5430
Immunogen	Synthetic peptide corresponding to residues surrounding S7 of human POLR2A CTD repeat YSPTSPS
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, ELISA, CHIP
Cross Reactivity	Human, Mouse, Rat
SWISS	P24928
Protein Weight	217kDa
Shipping	Ice bag