


## Anti-Phospho-Ppp1R12A (T853) Rabbit Polyclonal Antibody For Wb And Elisa - O14974

Item Number: CM0001046



### Introduction

Rabbit polyclonal antibody specific to phosphorylated T853 of PPP1R12A (MYPT1). Validated for WB and ELISA, reacts with human, mouse, and rat. Ideal for myosin phosphatase and cytoskeleton research.

[Learn More](#)

| Parameter        | Value                                                                          |
|------------------|--------------------------------------------------------------------------------|
| Product Name     | Phospho-PPP1R12A-T853 Rabbit pAb                                               |
| Remarks/Alias    | MBS; GUBS; M130; MYPT1; Phospho-PPP1R12A-T853                                  |
| Species          | Human                                                                          |
| GeneID (Human)   | 4659                                                                           |
| GeneID           | 4659                                                                           |
| Immunogen        | A synthetic phosphorylated peptide around T853 of human PPP1R12A (NP_002471.1) |
| Source           | Rabbit                                                                         |
| Category         | Phosphorylated Antibodies                                                      |
| Application      | WB, ELISA                                                                      |
| Cross Reactivity | Human, Mouse, Rat                                                              |
| SWISS            | O14974                                                                         |
| Protein Weight   | 115kDa                                                                         |
| Shipping         | Ice bag                                                                        |