

Anti-Phospho-Pten-S380 Rabbit Monoclonal Antibody For Wb, Elisa - P60484

Item Number: CM0000956



Introduction

Phospho-PTEN-S380 Rabbit mAb (P60484) validated for Western blot and ELISA; cross-reacts with mouse PTEN. Enables study of PTEN phosphorylation in PI3K-AKT signaling, tumor suppression, and cell migration.

[Learn More](#)

Parameter	Value
Product Name	Phospho-PTEN-S380 Rabbit mAb
Remarks/Alias	BZS; DEC; CWS1; GLM2; MHAM; TEP1; MMAC1; PTEN1; 10q23del; PTENbeta; Phospho-PTEN-S380
Species	Human
Gene ID	5728
Immunogen	Synthetic Peptide: A synthetic phosphorylated peptide around S380 of human PTEN (P60484)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	P60484
Protein Weight	47kDa
Shipping	Ice bag