


Anti-Phospho-Rb-S780 Rabbit Monoclonal Antibody For Ihc-P, Elisa - P06400

Item Number: CM0001144



Introduction

Phospho-specific rabbit monoclonal antibody targeting human RB1 phosphorylated at Ser780. Suitable for IHC-P and ELISA applications; reacts with human, mouse, and rat. Used in cell cycle and tumor suppressor research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-RB-S780 Rabbit mAb
Remarks/Alias	RB; pRb; OSRC; pp110; p105-Rb; PPP1R130; p110-RB1; Phospho-RB-S780
Species	Human
Gene ID (Human)	5925
Gene ID	5925
Immunogen	Synthetic Peptide: A synthetic phosphorylated peptide around S780 of human RB (P06400)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P06400
Protein Weight	106kDa
Shipping	Ice bag