


Anti-Phospho-Rb (Ser780) Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P06400

Item Number: CM0001041



Introduction

Rabbit polyclonal antibody targeting phosphorylated RB (Ser780). Validated in WB, IHC-P, ELISA. Reacts with human, mouse, rat RB. Suitable for cell cycle and tumor suppressor research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-RB-S780 Rabbit pAb
Remarks/Alias	RB; pRb; OSRC; pp110; p105-Rb; PPP1R130; p110-RB1; Phospho-RB-S780
Species	Human
GeneID (Human)	5925
GeneID	5925
Immunogen	Synthetic peptide
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P06400
Protein Weight	106kDa
Shipping	Ice bag