

Anti-Phospho-Ret (Tyr905) Rabbit Polyclonal Antibody For Wb & Elisa - P07949

Item Number: CM0000856



Introduction

Phospho-RET (Tyr905) rabbit polyclonal antibody (WB, ELISA) detects human RET phosphorylated at tyrosine 905, a receptor tyrosine kinase critical for cell proliferation, migration, and survival. Suitable for oncology and developmental biology studies.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Ret-Y905 Rabbit pAb
Remarks/Alias	PTC; MTC1; HSCR1; MEN2A; MEN2B; CDHF12; CDHR16; RET-ELE1; Phospho-Ret-Y905
Species	Human
GeneID (Human)	5979
GeneID	5979
Immunogen	Synthetic peptide / A synthetic phosphorylated peptide around Y905 of human Ret (NP_066124.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P07949
Protein Weight	124kDa
Shipping	Ice bag