

## Anti-Phospho-Shp1-Y564 (Ptpn6) Polyclonal Antibody For Wb, Elisa - P29350

Item Number: CM0001618



### Introduction

Rabbit polyclonal antibody targeting phosphorylated SHP-1 at tyrosine 564, validated for western blotting and ELISA in human and mouse; suitable for immune signaling and cancer research.

[Learn More](#)

| Parameter        | Value                                                                                           |
|------------------|-------------------------------------------------------------------------------------------------|
| Product Name     | Phospho-SHP1-Y564 Rabbit pAb                                                                    |
| Remarks/Alias    | HCP; HCPH; SHP1; SHP-1; HPTP1C; PTP-1C; SHP-1L; SH-PTP1; Phospho-SHP1-Y564                      |
| Species          | Human                                                                                           |
| GeneID (Human)   | 5777                                                                                            |
| GeneID           | 5777                                                                                            |
| Immunogen        | Synthetic peptide   A synthetic phosphorylated peptide around Y564 of human PTPN6 (NP_002822.2) |
| Source           | Rabbit                                                                                          |
| Category         | Phosphorylated Antibodies                                                                       |
| Application      | WB, ELISA                                                                                       |
| Cross Reactivity | Human, Mouse                                                                                    |
| SWISS            | P29350                                                                                          |
| Protein Weight   | 68kDa                                                                                           |
| Shipping         | Ice bag                                                                                         |