

Anti-Phospho-Srf (Ser103) Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P11831

Item Number: CM0001024



Introduction

Rabbit polyclonal antibody specific to phosphorylated Ser103 of human SRF. Validated for Western blot, immunohistochemistry (paraffin), and ELISA. Cross-reacts with mouse and rat. Useful for studying SRF-mediated transcription and cytoskeletal regulation.

[Learn More](#)

Parameter	Value
Product Name	Phospho-SRF-S103 Rabbit pAb
Remarks/Alias	MCM1; Phospho-SRF-S103
Species	Human
Gene ID	6722
Immunogen	A synthetic phosphorylated peptide around S103 of human SRF (NP_003122.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P11831
Protein Weight	52kDa
Shipping	Ice bag