

Anti-Phospho-Stat1 (S727) Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P42224

Item Number: CM0001952



Introduction

Rabbit polyclonal antibody recognizing phospho-STAT1 (Ser727). Validated in WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat STAT1. Ideal for interferon signaling and immune response studies.

[Learn More](#)

Parameter	Value
Product Name	Phospho-STAT1-S727 Rabbit pAb
Remarks/Alias	CANDF7; IMD31A; IMD31B; IMD31C; ISGF-3; STAT91; Phospho-STAT1-S727
Species	Human
GeneID (Human)	6772
Immunogen	Synthetic peptide (A synthetic phosphorylated peptide around S727 of human STAT1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P42224
Protein Weight	87kDa
Shipping	Ice bag