

Anti-Phospho-Stk3/Mst2-T384 Rabbit Polyclonal Antibody For Wb, Elisa - Q13188

Item Number: CM0001335



Introduction

Rabbit polyclonal antibody specific for phosphorylated Thr384 of human STK3/MST2. Validated for WB and ELISA; cross-reacts with human, mouse, and rat. Suitable for studying apoptosis and Hippo pathway signaling.

[Learn More](#)

Parameter	Value
Product Name	Phospho-STK3/MST2-T384 Rabbit pAb
Remarks/Alias	KRS1; MST2; Phospho-STK3/MST2-T384
Species	Human
Gene ID (Human)	6788
Gene ID	6788
Immunogen	Synthetic Peptide: A synthetic phosphorylated peptide around T384 of human STK3/MST2 (NP_006272.2)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13188
Protein Weight	56kDa
Shipping	Ice bag