

Anti-Phox2B Rabbit Polyclonal Antibody For Wb, Elisa - Q99453

Item Number: CM0022225



Introduction

Anti-PHOX2B rabbit polyclonal antibody for WB and ELISA detects human and mouse PHOX2B, a transcription factor critical for noradrenergic neuron development. Ideal for neuroblastoma and brain research.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated]PHOX2B Rabbit pAb
Remarks/Alias	CCHS; PMX2B; NBLST2; NBPhox; PHOX2B
Species	Human
GeneID (human)	8929
GeneID	8929
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 260-314 of human PHOX2B (NP_003915.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q99453
Protein Weight	32kDa
Shipping	Ice bag