

Anti-Pi3 Kinase P85 Alpha/Beta Rabbit Polymab® For Wb, If-P, Ihc-P, Elisa - P27986/O00459

Item Number: CM0007528



Introduction

Validated for WB, IF-P, IHC-P, and ELISA. Reacts with human, mouse, rat. Targets PIK3R1/PIK3R2 (p85 alpha/beta). Optimal for studying PI3K/AKT signaling in cancer, metabolism, and immune response research.

[Learn More](#)

Parameter	Value
Product Name	PI3 Kinase p85 alpha/beta Rabbit PolymAb®
Remarks/Alias	PIK3R1/PIK3R2
Species	Human
Gene ID	5295 (PIK3R1) / 5296 (PIK3R2)
Immunogen	Synthetic peptide; A synthesized peptide derived from human PI3 Kinase p85 alpha/beta
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P27986/O00459
Protein Weight	42-54 kDa/81-84 kDa
Shipping	Ice bag