


Anti-Pias1+Pias2 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - O75925/O75928

Item Number: CM0012172



Introduction

Rabbit monoclonal antibody recognizing PIAS1 and PIAS2, suitable for WB, IHC-P, and ELISA. Detects endogenous levels in human, mouse, and rat. Ideal for investigating SUMOylation, STAT signaling, and chromatin regulation.

[Learn More](#)

Parameter	Value
Product Name	PIAS1+PIAS2 Rabbit mAb
Remarks/Alias	DDXBP1; GBP; GU/RH-II; ZMIZ3; PIAS1+PIAS2
Species	Human
GeneID (human)	8554/9063
GeneID	8554/9063
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 552-651 of human PIAS1+PIAS2 (O75925).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O75925/O75928
Protein Weight	76kDa
Shipping	Ice bag