


Anti-Pik3C2A Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - O00443

Item Number: CM0018803



Introduction

Rabbit polyclonal antibody against human PIK3C2A (PI3K-C2 α), recognizing a ~191 kDa protein. Validated for WB, IHC-P, IF/ICC, IP, ELISA in human, mouse, and rat samples. Suitable for studies of PI3K signaling, insulin secretion, and intracellular trafficking.

[Learn More](#)

Parameter	Value
Product Name	PIK3C2A Rabbit pAb
Remarks/Alias	CPK; OCKSD; PI3-K-C2A; PI3K-C2alpha; PI3K-C2-alpha; PI3-K-C2(ALPHA); PIK3C2A
Species	Human
GeneID (Human)	5286
GeneID	5286
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human PIK3C2A (NP_002636.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O00443
Protein Weight	191kDa
Shipping	Ice bag