


Anti-Pim2 Rabbit Monoclonal Antibody For Wb, Ip, Elisa - Q9P1W9

Item Number: CM0007397



Introduction

Rabbit monoclonal antibody targeting human PIM2 serine/threonine-protein kinase, a proto-oncogene regulating cell survival and proliferation. Suitable for WB, IP, ELISA applications. Detects a band of approximately 34 kDa.

[Learn More](#)

Parameter	Value
Product Name	PIM2 Rabbit mAb
Remarks/Alias	PIM2; serine/threonine-protein kinase pim-2
Species	Human
GeneID (Human)	11040
GeneID	11040
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 200-300 of human PIM2 (Q9P1W9).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human
SWISS	Q9P1W9
Protein Weight	34kDa
Shipping	Ice bag