


Anti-Pitpnm1 Polyclonal Antibody For Wb, Elisa - 000562

Item Number: CM0016929



Introduction

Rabbit polyclonal antibody targeting human PITPNM1. Validated for WB and ELISA, cross-reacts with human, mouse, and rat. PITPNM1 is a phosphatidylinositol transfer protein essential for membrane trafficking and signaling. For research use only.

[Learn More](#)

Parameter	Value
Product Name	PITPNM1 Rabbit pAb
Remarks/Alias	Rd9; NIR2; RDGB; DRES9; RDGB1; RDGBA; PITPNM; RDGBA1; PITPNM1
Species	Human
GeneID (Human)	9600
GeneID	9600
Immunogen	Recombinant protein / Recombinant fusion protein containing a sequence corresponding to amino acids 110-190 of human PITPNM1 (NP_004901.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O00562
Protein Weight	135kDa
Shipping	Ice bag