


Anti-Pka C-Alpha Rabbit Monoclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - P17612

Item Number: CM0011889



Introduction

Recombinant rabbit monoclonal antibody against PKA C-alpha (PRKACA), validated for WB, IF/ICC, IHC-P, ELISA. Reacts with human, mouse, rat. Suitable for studying cAMP signaling, cancer, and platelet function.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] PKA C-alpha Rabbit mAb
Remarks/Alias	CAFD1; PKACA; PPNAD4
Species	Human
Gene ID (Human)	5566
Gene ID	5566
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-351 of human PKA C-alpha (PRKACA) (NP_002721.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P17612
Protein Weight	41 kDa/40 kDa
Shipping	Ice bag