

Anti-Pkc Zeta Polyclonal Antibody For Wb, If/Icc, Elisa - Q05513

Item Number: CM0029470

Introduction

Rabbit polyclonal antibody targeting PKC zeta (nPKC-zeta), validated for WB, IF/ICC, ELISA in human, mouse, rat. Useful for investigating PI3K pathway, MAPK cascade, cell polarity, NF-kB activation, insulin signaling, and inflammation.

[Learn More](#)

Parameter	Value
Product Name	PKC zeta Rabbit pAb
Remarks/Alias	PKC2; PKC-ZETA; PKC zeta
Species	Human
Gene ID (Human)	5590
Gene ID	5590
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 476-547 of human PKC zeta (NP_002735.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q05513
Protein Weight	68kDa
Shipping	Ice bag