


Anti-Pkn1 Polyclonal Antibody For Wb, Elisa - Q16512

Item Number: CM0013156



Introduction

Rabbit polyclonal PKN1 antibody validated for WB and ELISA. Reacts with human and mouse PKN1, a serine/threonine kinase involved in cytoskeletal regulation, cell migration, and transcription. Ideal for cancer and cell signaling research.

[Learn More](#)

Parameter	Value
Product Name	PKN1 Rabbit pAb
Remarks/Alias	DBK; PKN; PAK1; PRK1; PAK-1; PRKCL1; PKN-ALPHA; PKN1
Species	Human
Gene ID	5585 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human PKN1 (NP_002732.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q16512
Protein Weight	104kDa
Shipping	Ice bag