


Anti-Pla2G1B Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P04054

Item Number: CM0026613



Introduction

Rabbit polyclonal antibody against human PLA2G1B (phospholipase A2) for Western blot, IHC-P, and ELISA. Cross-reacts with human and mouse. Target is involved in phospholipid metabolism, inflammation, and anti-helminth defense. Ideal for digestive enzyme research.

[Learn More](#)

Parameter	Value
Product Name	PLA2G1B Rabbit pAb
Remarks/Alias	PLA2; PLA2A; PPLA2; PLA2G1B
Species	Human
GeneID (Human)	5319
GeneID	5319
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 23-148 of human PLA2G1B (NP_000919.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P04054
Protein Weight	16kDa
Shipping	Ice bag