

Anti-Placental Alkaline Phosphatase (Plap) Polyclonal Antibody For Wb, Ip, Elisa - P05187

Item Number: CM0019462



Introduction

CamelBio supplies a rabbit polyclonal antibody targeting human placental alkaline phosphatase (PLAP), suitable for Western blot, immunoprecipitation, and ELISA applications in placental biology and cancer research.

[Learn More](#)

Parameter	Specification
Product Name	Placental alkaline phosphatase (PLAP) Rabbit pAb
Remarks/Alias	ALP; IAP; ALPI; PALP; PLAP; PLAP-1; Placental alkaline phosphatase (PLAP)
Species	Human
GeneID (Human)	250
GeneID	250
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 260-460 of human Placental alkaline phosphatase (PLAP) (NP_001623.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human
SWISS	P05187
Protein Weight	58kDa
Shipping	Ice bag