


Anti-Plc Gamma 1 (Plcg1) Rabbit Polyclonal Antibody For Wb, Elisa - P19174

Item Number: CM0025594



Introduction

Rabbit polyclonal antibody targeting human PLCG1 (phospholipase C gamma-1), validated for WB and ELISA with cross-reactivity to mouse and rat. Essential tool for studying intracellular signaling cascades and actin cytoskeleton reorganization.

[Learn More](#)

Parameter	Description
Product Name	PLC gamma 1 (PLCG1) Rabbit pAb
Remarks/Alias	PLC1; NCKAP3; PLC-II; PLC148; PLCgamma1; PLC gamma 1 (PLCG1)
Species	Human
GeneID (Human)	5335
GeneID	5335
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 700-800 of human PLC gamma 1 (PLCG1) (NP_002651.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P19174
Protein Weight	149kDa
Shipping	Ice bag