


Anti-Plgf Rabbit Monoclonal Antibody For Wb, If-P, Elisa - P49763

Item Number: CM0006330



Introduction

High-quality rabbit monoclonal antibody against human PLGF (Placenta growth factor), validated for WB, IF-P, ELISA. Cross-reacts with mouse and rat. Ideal for angiogenesis and cancer research.

[Learn More](#)

Parameter	Value
Product Name	PLGF Rabbit mAb
Remarks/Alias	PGFL; PIGF; PLGF; PIGF-2; D12S1900; SHGC-10760
Species	Human
GeneID (Human)	5228
GeneID	5228
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human PLGF (P49763).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P49763
Protein Weight	25 kDa
Shipping	Ice bag