


Anti-Plgf Rabbit Polyclonal Antibody For Wb And Elisa - P49763

Item Number: CM0018951



Introduction

High-quality rabbit polyclonal antibody against human PLGF (Placenta growth factor), validated for WB and ELISA. Recognizes human, mouse, and rat PLGF. Ideal for angiogenesis and tumor growth research.

[Learn More](#)

Parameter	Value
Product Name	PLGF Rabbit pAb
Remarks/Alias	PGFL; PIGF; PLGF; PIGF-2; D12S1900; SHGC-10760
Species	Human
GeneID (Human)	5228
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 19-170 of human PLGF (NP_002623.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P49763
Protein Weight	25 kDa
Shipping	Ice bag