


Anti-Plk1 Rabbit Monoclonal Antibody (Plk1/Stpk13) For Wb, Elisa - P53350

Item Number: CM0003640



Introduction

High-quality rabbit monoclonal antibody against human PLK1, validated for Western blot and ELISA. Reacts with human, mouse, rat, and monkey PLK1. Ideal for cell cycle and mitosis research.

[Learn More](#)

Parameter	Value
Product Name	PLK1 Rabbit mAb
Remarks/Alias	PLK; STPK13; PLK1
Species	Human
Gene ID	5347
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 272-495 of human PLK1 (NP_005021.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat, Monkey
SWISS	P53350
Protein Weight	68kDa
Shipping	Ice bag