


Anti-Plod3 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - O60568

Item Number: CM0017036



Introduction

Rabbit polyclonal antibody against human PLOD3 (LH3) validated for WB, IHC-P, IF/ICC, and ELISA. Recognizes human and mouse PLOD3, a multifunctional procollagen lysine hydroxylase and glycosyltransferase essential for collagen biosynthesis and basement membrane formation.

[Learn More](#)

Parameter	Value
Product Name	PLOD3 Rabbit pAb
Remarks/Alias	LH3; PLOD3
Species	Human
GeneID (Human)	8985
GeneID	8985
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-280 of human PLOD3 (NP_001075.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	O60568
Protein Weight	85kDa
Shipping	Ice bag