

Anti-Pmepa1 Rabbit Polyclonal Antibody For Wb And Elisa - Q969W9

Item Number: CM0015462

Introduction

Rabbit polyclonal antibody targeting human PMEPA1 (TMEPAI), validated for WB and ELISA, cross-reactive with mouse and rat. Suitable for studying TGF-beta signaling and androgen receptor regulation in cancer research.

[Learn More](#)

Parameter	Value
Product Name	PMEPA1 Rabbit pAb
Remarks/Alias	STAG1; TMEPAI; PMEPA1
Species	Human
Gene ID	56937 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 173-252 of human PMEPA1 (NP_954638.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q969W9
Protein Weight	32kDa
Shipping	Ice bag