

Anti-Pnlip Rabbit Polyclonal Antibody For Wb, If-P, Elisa - P16233

Item Number: CM0030680



Introduction

Rabbit polyclonal antibody targeting human pancreatic lipase (PNLIP), validated for WB, IF-P, ELISA with cross-reactivity to mouse and rat. Ideal for lipid metabolism and pancreatic function research.

[Learn More](#)

Parameter	Specification
Product Name	PNLIP Rabbit pAb
Remarks/Alias	PL; PTL; PNLIPD; PNLIP
Species	Human
GeneID (Human)	5406
GeneID	5406
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 166-465 of human PNLIP (NP_000927.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P16233
Protein Weight	51kDa
Shipping	Ice bag