

Anti-Pnlip Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P16233

Item Number: CM0026654



Introduction

Rabbit polyclonal antibody targeting human pancreatic lipase (PNLIP). Validated for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with mouse. Ideal for studying fat digestion and metabolic disorders.

[Learn More](#)

Parameter	Value
Product Name	PNLIP Rabbit pAb
Remarks/Alias	PL; PTL; PNLIPD; PNLIP
Species	Human
Gene ID (Human)	5406
Gene ID	5406
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 166-465 of human PNLIP (NP_000927.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P16233
Protein Weight	51kDa
Shipping	Ice bag