


## Anti-Polr2A Rabbit Polyclonal Antibody For Wb, If/Icc, Ip, Elisa - P24928

Item Number: CM0024670



### Introduction

High-quality rabbit polyclonal antibody targeting human POLR2A (RPB1), the catalytic subunit of RNA polymerase II. Validated for WB, IF/ICC, IP, and ELISA. Cross-reacts with mouse and rat. Ideal for transcription mechanism research.

[Learn More](#)

| Parameter        | Value                                                                                                             |
|------------------|-------------------------------------------------------------------------------------------------------------------|
| Product Name     | POLR2A Rabbit pAb                                                                                                 |
| Remarks/Alias    | RPB1; RPO2; POLR2; POLRA; RPBh1; RPOL2; NEDHIB; RplILS; hsRPB1; hRPB220; POLR2A                                   |
| Species          | Human                                                                                                             |
| GeneID (Human)   | 5430                                                                                                              |
| GeneID           | 5430                                                                                                              |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human POLR2A (NP_000928.1) |
| Source           | Rabbit                                                                                                            |
| Category         | Polyclonal Antibodies                                                                                             |
| Application      | WB, IF/ICC, IP, ELISA                                                                                             |
| Cross Reactivity | Human, Mouse, Rat                                                                                                 |
| SWISS            | P24928                                                                                                            |
| Protein Weight   | 217kDa                                                                                                            |
| Shipping         | Ice bag                                                                                                           |